













Convenience and Security Options

		
<p>Wallmount 134201 -TAT-5000 not included</p>	<p>Quick Release Security System: 124277 (Core) 124278 (Oral Equiv) Includes TAT-5000, Wallmount 134201, Latex free coiled cable</p>	<p>Quick Release Security System: 124286 (Core) 124287 (Oral Equiv) Includes TAT-5000, Wallmount 134201, 8' Latex free vinyl covered steel cable</p>
		
<p>Dual Security System 124280 (Core) 124281 (Oral Equiv) Includes: TAT-5000, Requires Keyless Locking Wallmount 134307</p>	<p>Quick Release Security System: 124288 (Core) 124289 (Oral Equiv) Includes TAT-5000, Wallmount 134305, Latex free coiled cable</p>	<p>Quick Release Security System: 124290 (Core) 124291 (Oral Equiv) Includes TAT-5000, Wallmount 134305, 8' Latex free vinyl covered steel cable</p>
		
<p>Dual Security System 124280 (Core) 124281 (Oral Equiv) Includes: TAT-5000, Requires Keyless Locking Wallmount 134308</p>	<p>Quick Release Security System: 124292 (Core) 124293 (Oral Equiv) Includes TAT-5000, Wallmount 134306, Latex free coiled cable</p>	<p>Quick Release Security System: 124294 (Core) 124295 (Oral Equiv) Includes TAT-5000, Wallmount 134306, 8' Latex free vinyl covered steel cable</p>
 <p>Roll-A-Bout Part # 129455 5-leg rolling stand, double basket, 6' vinyl covered steel cable, lock and</p>	 <p>8' Security Cable</p> <p>Part #124309 - Coiled Cable Part #124311 - Latex Free Coiled Cable</p>	 <p>Vinyl Covered Steel Cable Latex Free</p> <p>6' Cable - Part #134302 8' Cable - Part #134030</p>

- **Quick Release System:** A discrete method for easily removal and reattachment of the cable when required.
- **Dual Security System:** A system requiring return of the instrument to the holder after 5 temperature measurements.
- **Suggestions for Selection:**
 - If instruments are to be secured at the nurses' station, consider the keyless locking wallmount or Dual Security System.
 - If instruments are to be mounted one per bed, or in a dedicated area (e.g. triage), consider Quick Release System.
 - If instruments are to be attached to an existing rolling stand, or a vital signs monitor, consider the 6' vinyl covered steel cable 134302.